

H5N1 NA / Neuraminidase (H274Y)

Catalog Number: 11676-VNAHC1



Sino Biological Inc.
Biological Solution Specialist

General Information

Gene Name Synonym:

NA

Protein Construction:

A DNA sequence encoding the Influenza A virus (A/Anhui/1/2005 (H5N1)) neuraminidase (ABU94738.1) (Met 1-Lys 449) was expressed, the cell lysates are collected, and bio-activity was tested

Source: Influenza A H5N1

Expression Host: Human Cells

QC Testing

Bio-activity:

Measured by its ability to cleave a fluorogenic substrate, 2'-(4-Methylumbelliferyl)- α -D-N-acetylneuraminic acid

The specific activity is > 300 U

One unit is defined as the amount of enzyme required to cleave 1 nmole of 2'-(4-Methylumbelliferyl)- α -D-N-acetylneuraminic acid per minute at pH 7.5 at 37°C.

Stability:

Samples are stable for up to twelve months from date of receipt at -70 °C

Molecular Mass:

The recombinant influenza A H5N1 Neuraminidase comprises 450 amino acids and has a predicted molecular mass of 49.2 kDa

Formulation:

Lyophilized from PBS, 0.6% Triton X-100, 7% Trehalose, 6% Mannitol, pH7.4

Please [contact us](#) for any concerns or special requirements.

Usage Guide

Storage:

Store it under sterile conditions at -70°C upon receiving. Recommend to aliquot into smaller quantities for optimal storage.

Avoid repeated freeze-thaw cycles.

Reconstitution:

It is recommended that 1 ml sterile water be added to the vial to prepare a stock solution.

Protein Description

Neuraminidase (NA) is a major membrane glycoproteins found on the surface of influenza virus. NA specifically catalyzes the hydrolysis removal of terminal sialic acid residues from viral and cellular glycoconjugates. It is known that HA binds to the sialic acid-containing receptors on the surface of host cells during initial infection, and at the end of an infectious cycle, NA cleaves the HA-sialic acid bondage from the newly formed virions and the host cell receptors during budding. NA thus is described as a receptor-destroying enzyme which facilitates virus release and efficient spread of the progeny virus from cell to cell. NA is a single-pass type II membrane protein which exists as a homotetramer, and the transmembrane domain is involved in lipid raft association during intracellular transport. NA is suggested to play a role in the determination of host range restriction on replication and virulence. Nine subtypes of NA have been identified, and subtypes N1 and N2 have been positively linked to epidemics in man.

Influenza A virus subtype H5N1, also known as "bird flu", A(H5N1) or simply H5N1, is a subtype of the Influenza A virus which can cause illness in humans and many other animal species. H5N1 is easily transmissible between birds facilitating a potential global spread of H5N1. It is mainly spread by domestic poultry, both through the movements of infected birds and poultry products and through the use of infected poultry manure as fertilizer or feed. Humans with H5N1 have typically caught it from chickens, which were in turn infected by other poultry or waterfowl.

References

1. Barman, S. et al., 2000, J. Virol. 74: 6538-45.
2. Colman, P.M. et al., 1983, Nature. 303: 41-4.
3. Suzuki, T. et al., 2005, J. Virol. 79: 11705-15.
4. Shinya K, et al., 2006, Nature. 440 (7083): 435-6.
5. von, Itzstein, M. 2007, Nat. Rev. Drug. Discov. 6: 967-74.
6. Christophe F, et al., 2009, Science. 324:1557-61.

Manufactured By Sino Biological Inc., FOR RESEARCH USE ONLY. NOT FOR USE IN HUMANS.

Fax :+86-10-51029969 ● Tel:+86-400-890-9989 ● <http://www.sinobiological.com>